# Silagum®-Comfort

908170/08.2008

## Gebrauchsinformation

## Deutsch

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## Instructions for use

English

Permanently soft VPS-based relining material

### Indication

Relining material for total and partial dentures

### Recommended use

#### A. Direct relining procedure in the mouth

##### 1. Preparing the denture

If necessary, remove old permanently soft relining material. Trim off a minimum of 1mm of all the areas needing to be relined. A rounded preparation of the denture ridges is recommended. Following this preparation protect the denture surface from moisture, saliva, and monomer containing resins.

##### 2. Applying the primer

Apply an even layer of primer and let the solvent evaporate completely (approx. 1 minute). Be sure to fully cover the marginal ridges with primer. Apply at least two coats. Firmly close the bottle after use to prevent evaporation of the solvent.

##### 3. Applying relining material

First apply the relining material on the ridges and then coat the base of the denture. The material layer should be at least 2mm thick. The total working time of Silagum-Comfort is 1:45 minutes at room temperature.

##### 4. Placing the denture into the mouth

Control moisture. Insert the denture into the patient's mouth no later than 1 minute after the Silagum-Comfort application and let the patient close in normal occlusion. Let the patient carry out gentle swallowing and masticatory movements (without exerting excessive pressure) with the mouth closed for 5 minutes.

Be sure to prevent uncured silicone from entering the patient's throat and esophagus. The functional movements should therefore be carried out with the patient sitting upright or slightly bent forward.

##### 5. Removing the denture from the mouth

The denture can be removed from the mouth after approx. 6 minutes. Clean and disinfect the denture. Wait an additional 10 minutes to trim and polish.

The dry and fat free denture can be corrected immediately by applying and modeling additional material. Let the denture set in the mouth for 7 minutes or 5 minutes in a pressure pot (40–50 °C).

##### 6. Finishing and varnishing

Remove excess material with sharp scissors or a scalpel and smooth the edges with a suitable silicone trimmer.

Apply the enclosed varnish to obtain a perfectly smooth surface. Mix the varnish evenly 1:1, apply bubble-free with a brush, and let dry for 15 minutes.

#### B. Indirect relining at the laboratory

##### 1. Pouring and mounting the denture

Following a functional impression the denture is poured with gypsum toward the jaw and mounted on a relining jig.

Do not use polyether impression materials! The impression material used should feature extended working time and appropriate viscosity. The application of a suitable separating medium (gypsum against resin) is recommended.

##### 2. Fabricating a mother matrix mold

Produce a negative by pouring the occlusal part of the denture or rows of teeth, which must be firmly fixed in its vertical relationship to the model. For this purpose, the use of a special relining device ensuring exact reduction of model and negative is recommended.

##### 3. Preparing the denture

After all the impression material has been removed trim the denture surface toward the jaw at least 1mm and clean. Protect the denture surface from moisture, saliva, and monomer containing resins.

##### 4. Applying the primer

Apply an even layer of primer and let the solvent evaporate completely (approx. 1 minute). Be sure to completely cover the marginal ridges. Apply at least two coats. Firmly close the bottle after use to prevent evaporation of the solvent.

##### 5. Applying relining material

First, apply the relining material on the ridges and then coat the base of the denture. The material layer should be at least 2mm thick. The total working time of Silagum-Comfort is 1:45 minutes at room temperature.

##### 6. Curing the denture on the model

Reposition the denture on the model, and, using the device, press the negative to the stop, and let set for 30 minutes (or 10 minutes in the pressure pot at 40–50 °C).

The dry and fat free relining can be corrected immediately by applying additional material and modeling it. Let the denture cure now for 7 minutes in the mouth or 5 minutes in a pressure pot (40–50 °C).

##### 7. Finishing and varnishing

Remove excess material with sharp scissors or a scalpel and smooth the edges with a suitable silicone trimmer. Apply the supplied varnish to obtain a perfectly smooth surface. Mix the varnish evenly 1:1, and apply with a brush – be sure to avoid air entrapments – and let dry for 15 minutes.

### Please Note

Do not expose relining to mechanical stress for 24 hours (e. g. cleaning with a brush).

Avoid contact of the relining material with latex gloves. The use of PVC or PE gloves is recommended.

### Specifications

	Silagum-Comfort	Silagum-Comfort Lack	Silagum-Comfort Primer
Working time*	1:45 min	–	2 min
Setting time in mouth	6 min	–	–
Setting time*	15 min	–	15 min
Evaporation time	–	1 min	–

\*Times refer to normal room temperature (23 °C/73 °F) and normal relative humidity (50%) conditions. Higher temperatures shorten, lower temperatures increase these times.

### Disinfection

The Automix-Dispenser and the trimmer and brush holders are reusable following disinfection but should be replaced in case of visible wear (e.g. scratches). Use a commercially available wipe or dipping solution. We recommend using only RKI (Robert Koch Institut) listed disinfection solutions.

The Automix-Dispenser and brush holder are also autoclavable.

### Composition

Silagum-Comfort: Addition-cured vinylpolysiloxanes, hydrogenpolysiloxanes, fillers, pigments, additives, platinum catalyst

Silagum-Comfort Primer: Ethylacetate, modified polyacrylate, additives

Silagum-Comfort Lack: Vinyl silicone, hydrogen silicone, additives, platinum catalyst

### Storage

Store and use in dry location at normal room temperature not exceeding 25 °C/77 °F. Do not store the material in the refrigerator or freezer. Leave used tips on the cartridge after use. They serve as a cap.

Do not use after expiration date.

### Packaging

#### Introductory package

1 Cartridge @ 25ml paste, 5 Automix-Tips, 1 Automix-Dispenser  
1 Bottle primer @ 5ml, 2 Bottles varnish @ 10ml  
25 Brushes, 5 Trimmers, 1 Holder REF 909087

<b>Refill</b>	
1 Cartridge @ 25ml paste, 5 Automix-Tips,	REF 909088
1 Cartridge @ 50 ml paste, 12 Automix-Tips,	REF 909951
<b>Value pack</b>	
8 Cartridges @ 50ml paste	REF 909950
<b>Accessories</b>	
2 Bottles @ 10ml Silagum-Comfort Lack	REF 909089
1 Bottle @ 5ml Silagum-Comfort Primer	REF 909090
1 Automix-Dispenser type 25 1:1	REF 110253
1 Automix-Dispenser type 50 1:1	REF 909507
50 Automix-Tips, Yellow 1:1	REF 909200
Federal Law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.	
Keep away from children! For dental use only!	

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